CEPHA ONE

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Telefone: +55 11 40231749 | e-mail: contact@biocontinental.com

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1- International Committee For Standardization in Hematology Thrombosis and Haemostasis. Amer. J. Clin. Path.

2- National Committee For Clinical Laboratory Standards: Colection, Transport And Processing Of Blood Specimens For Coagulation Test And Performance Of Coagulation Assays. Approved Guideline Fourth Edition. NCCLS document H12-A4, 2003

3- Banez,el, D.a., Koepek, Laboratory Monitoring Of Heparin Therapy. Effect Of Different Salts Of Heparin On The Activated Partial Thromboplastin Time. Amer. J. Clin. Path.

4- ADCOCK, D. M.; KRESSIN, D. C. Effect of 3.2% vs. 3.8% sodium citrate concentration on routine coagulation testing. AM. J. Clin. Pathol. 107(1):105-10, 1997.

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CEPHA ONE

CLINICAL SIGNIFICANCE

The Activated Partial Thromboplastin Time (APTT) is a test that monitors both the intrinsic pathway and a common pathway of coagulation and is useful in monitoring of Heparin therapy and to verify the deficiency all coagulation factors except factor VII. Are widely used in deficiencies identified in stage 1 of the coagulation factors VIII, IX, XI, XII, Fletcher factor.

PURPOSE AND BASIS OF THE METHOD

Method for measuring manual or automatic time of Activated Partial Thromboplastin Time in citrated plasma. For in vitro diagnostics use only. Consist on the determination of time of citrated plasma after adding reagents that contain a Plasmatic Activator (Ellagic Acid) and Phospholipids, that will act as substitutes to platelets. Initially plasma is incubated with the Contact Activator, at 37°C for 3 minutes. Following, add the Clot Forming (Calcium Chlorite) for recalcification, triggering the intrinsic coagulation mechanism. This step should be timed until the formation of the clot.

REAGENTS PROVIDED

Number 1 - Contact Activator (Activated Cephalin) - Store between 2 and 8°C. Contains: Rabbit brain extract, Ellagic Acid and preservative.

Number 2 - Clot Former (Calcium Chlorite Solution) - Store between 2 and 8°C. Contains: Calcium Chloride and preservative.

Presentation ready for use.

MATERIALS REQUIRED AND NOT PROVIDED

Coagulometer, 37°C water bath, clock or stopwatch, pipettes, test tubes, control plasmas. These items are found in the specialized market for articles for Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

Storage temperature should be between 2 and 8oC. Transportation under temperatures between 15 and 30oC should exceed 72 (seventy two) hours. Protect from light and avoid moisture. Do not freeze.

PRECAUTIONS

- 1- For professional in vitro diagnostic use only.
- 2- Strictly follow the proposed methodology to obtain exact results.
- 3- The water level in the water bath must be higher than the level of the reagents in the test tubes.
- 4- Avoid heating the reagent that will not be used. Reagent that is already heated should not be returned to the original bottle.
- 5- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.
- 6- To obtain information related to biosafety or in case of accidents with the product, consult the Product Safety Data Sheet.
- 7- Do not use the product in case of damage to the packaging.
- 8- It is essential that the equipment used is properly calibrated and subjected to periodic maintenance.

SAMPLES

Material: Plasma

Anticoagulant: Use Sodium Citrate 3.8% (0.130 M) or 3.2% (0.109 M) as anticoagulant, in the proportion of 9 parts of Blood to 1 of Anticoagulant (as an example: 4.5 mL of Blood + 0.5 mL of Anticoagulant), or commercial Citrate kits, according to the manufacturer's Instructions for Use, or vacuum tubes with citrate, available on the market.

Collection: Blood must be obtained through veinpuncture, avoiding hemolysis, prolonged tourniquet, bubble formation and tissue liquid aspiration. When evaluating only the coagulation test, always collect 2 tubes of citrate and discard the first one.

Preparation: Blood should be mixed with anticoagulant right after collection. Centrifuge immediately at 3000 rpm, during 15 minutes.

Remove plasma without pipetting red cells and the yellow layer.

Samples must be tested in less than 3 hours. If tests can not be made during this period plasma should be frozen for at maximum 1 week at -20°C.

POST OPEN STABILITY

The calibration stability of the CEPHA ONE kit after opening is 14 days. This stability may vary according to the test, equipment and environment conditions. Therefore, it is suggested to monitor the product's performance using control plasmas.

PROCEDURE

- 1- Preheat the Reagent N° 2 to 37°C. This procedure should be performed only with the volume that will be used to perform the tests. The reagent must remain warm for the entire procedure.
- 2- Separate the reaction tubes that will be used to perform the tests.
- 3- Pipette 100 µL of the Sample/Control into the respective reaction tube and preheat to 37°C for 2 minutes in a water bath.

4- Pipette 100 µL of Reagent N° 1 into the Sample/Control tube and incubate at 37°C, for 2 to 3 minutes.

5- Pipette 100 µL of Reagent N° 2 (preheated to 37°C) in the Sample/Control tube, simultaneously triggering the timer.

6- Leave the tube in the water bath at 37°C, shaking gently. Keep in the bath for about 20 seconds.

7- Remove the tube from the bath, gently tilt it once a second and stop the timer at the time of the clot formation. Note the clotting time.

CALCULATIONS

Calculate the average between coagulation time of each samples/controls.

Release all results in seconds.

CURVE OF SENSIBILITY TO HEPARIN

This procedure is useful to monitor therapies with Heparin.

A curve of sensibility to Heparin must be built adding the quantity of Heparin known to a mixture of normal plasmas (pool).

As a suggestion the following values of Heparin can be used:

0.0 / 0.1 / 0.2 / 0.4 / 0.6

Determine the dilution time dosages, plotting in a graphic the times and concentrations.

Each laboratory must establish its own calibration curve.

PROCEDURE LIMITATIONS

The concentration of citrate may interfere with coagulation time because it directly affects the concentration of calcium. Thus, APTT values are higher in tubes of 3.8% than in tubes of 3.2%⁴. Citrate at 3.8% concentration will release results higher than 3.2%. Therefore, it is necessary that the laboratory standardize so that the monitoring of the patient's diagnosis is always performed with the same anticoagulant concentration.

The concentration of sodium citrate should be adjusted in patients with hematocrit values above 55%. In samples with high hematocrit, the blood/anticoagulant ratio (9:1) is not maintained and may cause excess citrate to the volume of blood present in the tube. To calculate the citrate volume, use the formula described in item 4 of the Laboratory Diagnostic Manual for Hereditary Coagulopathies and Plaquetopathies of the Ministry of Health⁵.

REFERENCE VALUES

Interval should be used only as guidance. It is recommended that each laboratory establish, in population served, its own reference range, using a pool of normal plasmas or reference plasmas.

Age	Interval (seconds)
2 months	24 – 47
5 months	24 – 46
Children older than 6 months and Adults	24 – 38

The results provided by this CEPHA ONE kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

INTERFERENCES

Icteric, lipemic or hemolyzed samples.

Usage of plasma containing Heparin, EDTA or Sodium Oxalate.

Conditions and substances that may alter TTPA: Oral contraceptives, Estrogen, pregnancy, drugs Coumarin, Heparin, Asparaginase, Acetylsalicylic Acid, Metronidazole, Digitalis and Tetracycline.

PRODUCT PERFORMANCE

COMPARISON OF METHODS

The CEPHA ONE kit was compared with another commercially available kit for determining the Partial Activated Thromboplastin Time. 25 analyzes were performed and the results were evaluated. The linear equation obtained was $y = 1.0633x + 1.7026$ and the correlation coefficient 0.9891. With these results, it can be concluded that the kit has good methodological specificity.

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 2 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Time (sec)	31.330	39.680
Standard Deviation (sec)	0.359	1.026
Coefficient of Variation (%)	1.15	2.59

REPRODUCIBILITY

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 1 sample, obtaining the following results:

	Sample 1
Time (sec)	39.763
Standard Deviation (sec)	0.725
Coefficient of Variation (%)	1.82

SENSITIVITY

The sensitivity of CEPHA ONE kit to factors deficiency has been evaluated. The evaluation has shown that the kit is sensitive to factors VIII, IX, XI, XII e Prekallikrein.

SPECIFICITY

The specificity of the kit is guaranteed by the addition of one single activation source of coagulant factors, Ellagic Acid.

PRESENTATION

Presentation	Reagent 1	Reagent 2
1	6 x 2,5 mL	4 x 4 mL
2	6 x 5 mL	8 x 4 mL
3	12 x 2,5 mL	8 x 4 mL
4	12 x 5 mL	16 x 4 mL
5	6 x 5 mL	6 x 5 mL
6	6 x 5 mL	10 x 10 mL
7	10 x 5 mL	10 x 5 mL
8	10 x 5 mL	10 x 10 mL
9	10 x 10 mL	10 x 10 mL

CONSUMER INFORMATION:

Biocontinental guarantees the quality of the product, provided that all conditions are respected, and the storage and transport care, specified in the instruction for use, until the expiration date indicated on the packaging.

If further information or guidance is required, the customer should contact Biocontinental Biological Products.

Technical Assistance and request for Instruction of Use in printed format:

Phone: +55 11 40231749 | E-mail: contact@biocontinental.com

BIBLIOGRAPHIC REFERENCES

- 1- International Committee For Standardization in Hematology Trombosis and Haemostasis. Amer. J. Clin. Path.
- 2- National Committee For Clinical Laboratory Standards: Colection, Transport And Processing Of Blood Specimens For Coagulation Test And Performance Of Coagulation Assays. Approved Guideline Fourth Edition. NCCLS document H12-A4, 2003
- 3- Banez,el, D.a., Koepek, Laboratory Monitoring Of Heparin Therapy. Effect Of Different Salts of Heparin On The Activated Partial Thromboplastin Time. Amer. J. Clin. Path.
- 4- ADCOCK, D. M.; KRESSIN, D. C. Effect of 3.2% vs. 3.8% sodium citrate concentration on routine coagulation testing. AM. J. Clin. Pathol. 107(1): 105-10, 1997.
- 5- Manual de Diagnóstico Laboratorial das Coagulopatias Hereditárias e Plaquetopatias, Ministério da Saúde, 2010; Fisiologia da coagulação, anticoagulação e fibrinólise.
- 6- Rendrik F. Franco. Simpósio: Hemostasia e trombose. Capítulo I. 34: 229-237, jul./dez. 2001.

Manufacturer | Regularized by: Continental Produtos Biológicos Ltda
 Address: Rua Santana, 305, 13300-220 Itu- SP
 Site: www.biocontinental.com
 E-mail: contact@biocontinental.com
 CNPJ: 61.058.400/0001-42
 Brazilian Industry

ANVISA: 80132390003

Attention user: check if the version obtained from the Instruction for Use not printed corresponds to the product purchased.

Version: JUN/2025

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED