TT ONE

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TT ONE

CLINICAL SIGNIFICANCE

Thrombin Time determines the time that fibrinogen is converted to fibrin in the presence of thrombin. Abnormality of clotting time in TT can be caused when fibrinogen concentrations are below 70 to 100 mg / dL, dysfibrinogenemia, disseminated intravascular coagulation (DIC), liver disease, in some paraproteinemias. Clotting times are heavily influenced by heparin and hirudin, or any treatment or disease that interferes with the conversion of fibrinogen to fibrin.

PURPOSE AND BASIS OF THE METHOD

Method for manual or automated determination of Thrombin Time in citrated plasma. Only for *in vitro* diagnostic use.

It consists of measuring the plasma clotting time after the addition of thrombin (activated factor II), promoting the conversion of fibrinogen to fibrin.

REAGENTS PROVIDED

Number 1 - Thrombin Solution - Store at 2 to 10 °C. Contains bovine thrombin <6.0 NIH / mL, Calcium chloride (10 - 50 mmol / L), stabilizers and preservative.

Presentation ready for use.

MATERIALS REQUIRED AND NOT PROVIDED

Tips, pipettes, tubes, semi-automatic or automatic coagulation equipment (coagulometer), or timer (manual technique). These articles are available in the specialized market for Clinical Analysis Laboratories.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 10°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep out of the light. Do not freeze.

PRECAUTIONS

- 1- For professional in vitro diagnostic use only.
- 2- Strictly follow the proposed methodology to obtain exact results.
- 3- The water level in the water bath must be higher than the level of the reagents in the test tubes.
- 4- Avoid heating the reagent that will not be used. Reagent that is already heated should not be returned to the original bottle.
- 5- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.
- 6- To obtain information related to biosafety or in case of accidents with the product, consult the Product Safety Data Sheet.
- 7- Do not use the product in case of damage to the packaging.
- 8- It is essential that the equipment used is properly calibrated and subjected to periodic maintenance.

SAMPLES

Material: Plasma

Anticoagulant: Use Sodium Citrate 3.8% (0.130 M) or 3.2% (0.109 M) as anticoagulant, in the proportion of 9 parts of Blood to 1 of Anticoagulant (as an example: 4.5 mL of Blood + 0,5 mL of Anticoagulant), or commercial Citrate kits, according to the manufacturer's Instructions for Use, or vacuum tubes with citrate, available on the market.

Collection: Blood must be obtained by venipuncture, avoiding hemolysis, prolonged tourniquet, bubble formation and tissue fluid aspiration. When evaluating only the coagulation test, always collect 2 tubes of citrate and discard the first.

Preparation: The blood should be mixed with the anticoagulant immediately after collection.

Centrifuge immediately at 3000 rpm for 15 minutes.

Remove the plasma without pipetting red cells and the yellow layer.

The samples must be tested in less than 3 hours. If the test cannot be performed during this period, the plasma should be frozen for a maximum of 1 week at -20°C.

POST OPEN STABILITY

The stability of the TT ONE kit after opening is 14 days. This stability may vary according to the test conditions, the equipment and the environment. Therefore, it is suggested to monitor the product's performance using control plasmas.

PROCEDURE

Allow reagent N° 1 to reach room temperature before use.

1- Incubate 0.2 mL of plasma in a water bath at 37°C for 2 minutes.

2- Add 0.2 mL of unheated N° 1 reagent, simultaneously start the timer and record clot formation times.

RESULTS

The results of the Thrombin Time test are expressed in seconds. When the difference in the clotting times of the same sample is greater than 5%, it is recommended to repeat the test discarding the previous measurements.

TRACEABILITY

The TT ONE kit is traceable to the reference material NIBSC WHO 2nd International Standard for Thrombin 01/580 US FDA / CBER Thrombin Standard Lot K.

REFERENCE VALUES

Normal Values: <21 seconds.

The range should be used as a guide only. It is recommended that each Laboratory establish, in the population served, its own reference intervals

Aspects such as: quality of collection, conservation of the sample and equipment used can impact the sample's clotting time. It is recommended that each laboratory define its reference values for its population served.

INTERFERENT

According to the scientific literature, some substances are known as possible interferents for the coagulometry methodology for measuring thrombin time, which are causes of falsely elevated or decreased results.

Hemolysis and lipemia are interferents that, in excess, can change results unpredictably.

PRODUCT PERFORMANCE

COMPARISON OF METHODS

The TT ONE kit was compared with another method for evaluating commercially available Thrombin Time. 45 analyzes were performed and the results were evaluated.

The correlation coefficient equal to 0.9884. With these results, it can be concluded that the kit has good methodological specificity.

SPECIFICITY

The specificity of the kit is guaranteed by the addition of an exclusive source of thrombin, a specific activator of factor I (Fibrinogen), which eliminates the possibility of direct activation of other factors in the coagulation cascade.

REPEATABILITY

Repeatability was calculated from 10 consecutive determinations of three samples with different clotting times, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Time (sec)	18,48	22,32	19,87
Standard Deviation (sec)	0,17	0,20	0,26
Coefficient of Variation (%)	0,9	0,9	1,3

REPRODUCIBILITY

Reproducibility was calculated from 10 consecutive determinations of three samples with different clotting times, for 3 consecutive days, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Time (sec)	15,647	18,606	56,271
Standard Deviation (sec)	0,281	0,295	1,232
Coefficient of Variation (%)	1,80	1,58	2,19

SENSITIVITY

Kit sensitivity is defined by the concentration of thrombin coagulation factor activity in the kit, which is <6.0 NIH / mL.⁶

PRESENTATION

R1 - 2 x 2 mL

R1 - 5 x 2 mL

R1 - 6 x 2 mL

R1 - 10 x 2 mL

R1 - 12 x 2 mL

R1 - 6 x 5 mL

CONSUMER INFORMATION:

Biocontinental guarantees the quality of the product, provided that all conditions are respected, and the storage and transport care, specified in the instruction for use, until the expiration date indicated on the packaging.

If further information or guidance is required, the customer should contact Biocontinental Biological Products.

Technical Assistance and request for Instruction of Use in printed format:

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Brazilian Industry

ANVISA: 80132390002

Attention user: check if the version obtained from the Instruction for Use not printed corresponds to the product purchased.

Version: JUN/2025

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED