## FAMÍLIA DIMER ONE PLUS

### SIGNIFICADO DIAGNÓSTICO

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## DIMER ONE PLUS FAMILY

### DIAGNOSTIC SIGNIFICANCE

D-Dimer is a fibrin degradation product formed during activation of the coagulation system. Its measurement is used to rule out thromboembolic disease in outpatients with suspected deep vein thrombosis (DVT) and pulmonary embolism (PE). DVT and PE can cause sudden, fatal embolic events in the pulmonary arteries and other regions. Measuring plasma D-Dimer levels has been used as a screening strategy for subclinical DVT. ATVP is a high risk factor for stroke due to advanced age, hemiplegia, and coagulation disorders. Therefore, it is important to monitor the incidence level of D-dimer and the characteristics of VT in patients with acute stroke. D-dimer is also an important prognostic indicator of heart disease, playing an important role in monitoring the clinical picture after treatment and post-therapeutic evaluation of patients.

### PURPOSE AND METHOD

Method for the quantitative determination of D-dimer in a biological sample (plasmacitrated), using an immunoturbidimetric test. Only for in vitro diagnostic use. Latex particles, coated with anti-D-dimer antibodies, are agglutinated by the D-dimer present in the sample. The increase in turbidity caused by the agglutination of the latex particles is proportional to the concentration of D-dimer present in the sample and can be measured spectrophotometrically.

### REAGENTS PROVIDED

Reagent Number 1 - Buffer - Store between 2 and 8°C. Contains: Buffer <200 mmol/L, stabilizers, surfactant and preservative.

Reagent Number 2 - Latex D-Dimer - Store between 2 and 8°C. Contains: Suspension of latex particles sensitized with anti-D-dimer antibody, stabilizers and preservative.

Reagent Number 3 - Diluent Solution - Store between 2 and 8°C. Contains: Buffer <200 mmol/L, 0.9% NaCl and preservative.

Reagent Number 4 - Calibrator - Store between 2 and 8°C - Contains: D-dimer, stabilizers and preservative. Potentially infectious.

Reagents Number 5 and 6 – Controls Level 1 and 2 - Store between 2 and 8°C - Contains: D-dimer, stabilizers and preservative. Potentially infectious.  
Attention: The concentration of D-dimer in calibrators and controls varies according to the batch - As per bottle label.

### REQUIRED AND NOT PROVIDED MATERIALS

Automatic coagulation equipment or automatic biochemical equipment. These items are available on the market from suppliers specialized in Clinical Analysis Laboratories.

### STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be between 2 to 8°C. Transportation at temperatures up to 30°C should not exceed 5 days. Keep protected from light. Do not freeze..

### PRECAUTIONS

- 1- For professional in vitro diagnostic use only.
- 2- Follow the proposed methodology rigorously to obtain accurate results.
- 3- Avoid heating reagent that will not be used. Reagent that has already been heated should not be returned to the original bottle.
- 4- We recommend applying local, state, and federal environmental protection regulations so that the disposal of reagents and biological material is carried out in accordance with current legislation.
- 5- For information related to biosafety or in case of accidents with the product, check the Product Safety Data Sheet.
- 6- Do not use the product if the packaging is damaged.
- 7- It is essential that the equipment used is properly calibrated and subject to periodic maintenance

### SAMPLES

**Material:** Free citrate plasma with hemolysis and intense lipemia. The analyte is stable for 8 hours at 20-25°C, 4 days between 2 and 8°C, and 6 months at -20°C\*.

### STABILITY AFTER OPENING

The stability of DIMER ONE PLUS after opening is 30 days. The DIMER ONE PLUS calibrators and controls after opening and in use are stable for 5 days between 2 and 8°C and 20 days at -20°C, provided they are stored in a suitable container. This stability can vary according to the test conditions, the equipment, and the environment. Therefore, it is recommended to monitor the performance of the product using control samples.

### PROCEDURE

Ready-to-use reagents.  
Thoroughly homogenize Reagent No. 2 before use.

For reaction calibration, use the DIMER ONE PLUS Calibrator. For control, use the DIMER ONE PLUS Control.

The kit is intended for use only in automatic biochemical analyzers and/or automatic coagulation analyzers. Contact us to obtain the protocol for your automatic equipment.

SUGGESTED GENERAL PROGRAMMING PARAMETERS FOR AUTOMATIC BIOCHEMICAL EQUIPMENT	
Filter	570nm
Incubation Temperature	37°C – 98.6°F
Unit	µg/mL
Sample Volume R1 Volume* R2 Volume*	37,5µL 150µL 75 µL Incubation time may vary according to the equipment
Calibrator Calibration Mode	DIMER ONE PLUS Calibrator* Log- loglinear-(may vary according to the equipment)

\* The volumes of the reagents may vary, as long as they are proportional, according to the equipment.

GENERAL PARAMETERS OF SUGGESTED PROGRAMS FOR AUTOMATIC COAGULATION EQUIPMENT	
Filter	570nm
Incubation Temperature	37°C - 98.6°F
Unit	µg/mL
Sample Volume R1 Volume* R2 Volume*	25µL 100µL 50µL Incubation time may vary according to the equipment.
Calibrator Calibration Mode	DIMER ONE PLUS Calibrator* Log-loglinear(may vary according to the equipment)

\*The volumes of the reagents may vary, as long as they are proportional, according to the equipment.

**CALIBRATION:** Dilute the calibrator in the following proportions: 1/1, 1/2, 1/4, 1/8, and 1/16, using Reagent number 3 (Diluent Solution).

### RESULTS

DIMER ONE PLUS results are expressed in µg/mL.

### TRACEABILITY

The calibrators and controls of the kit were normalized using a purified preparation of D-dimer protein in mass and concentration.

### REFERENCE VALUES

Normal Values: < 0.5 µg/mL.  
The range should be used for guidance only. It is recommended that each laboratory establish its own reference ranges for the population served.

### INTERFERENTS

Normal Values: < 0.5 µg/mL.  
The range should be used for guidance only. It is recommended that each laboratory establish its own reference ranges for the population served..

### INTERFERENTS

According to the scientific literature, some substances are known as possible interferents for the coagulometry methodology used to measure thrombin time, which can cause falsely elevated or decreased results.

### DOSAGE FORM

#### DIMER ONE PLUS

R1: 1 x 10,0 mL / R2: 1 x 7,0 mL / R3: 1 x 10,0 mL  
R1: 1 x 7 mL / R2: 1 x 4 mL / R3: 1 x 7 mL  
R1: 1 x 5,0 mL / R2: 1 x 2,0 mL / R3: 1 x 6,0 mL

#### DIMER ONE PLUS + CALIBRATOR + CONTROL

R1: 1 x 10,0 mL / R2: 1 x 7,0 mL / R3: 1 x 10,0 mL / R4: 1 x 1 mL / R5: 1 x 1 mL / R6: 1 x 1 mL  
R1: 1 x 5,0 mL / R2: 1 x 2,0 mL / R3: 1 x 6,0 mL / R4: 1X1 mL / R5: 1x1 mL / R6: 1x1 mL  
R1: 2 x 10,0 mL / R2: 2 x 2,0 mL / R3: 2 x 8,0 mL / R4: 1X1 mL / R5: 2x1 mL / R6: 2x1 mL

### PRODUCT PERFORMANCE

#### COMPARISON OF METHODS

DIMER ONE PLUS kit It was compared with another commercially available method for evaluating D-dimer. Analyses were performed and the results were evaluated. The linear equation  $y = 1.0021x + 0.0339$  and the correlation coefficient 0.9995 were obtained. With these results, it can be concluded that the kit has good methodological specificity.

#### REPEATABILITY

Repeatability was calculated from 20 consecutive determinations of two samples with different concentrations, obtaining the following results.

Repeatability	Sample 1	Sample 2
Average concentration (µg/mL)	2,58	9,56
Standard Deviation (µg/mL)	0,04	0,07
Coefficient of Variation (%)	1,56%	0,69%

#### REPRODUCIBILITY

Reproducibility was calculated from 10 consecutive determinations of two samples with different concentrations over 3 consecutive days, obtaining the following results.

Reproducibility	Sample 1	Sample 2
Average concentration (µg/mL)	2,56	9,59
Standard Deviation (µg/mL)	0,04	0,07
Coefficient of Variation (%)	1,60%	0,71%

#### ANALYTICAL SENSITIVITY

Sensitivity was defined as the lowest concentration of D-dimer that presents a coefficient of variation less than or equal to 20%. Considering a coefficient of variation of 20%, the reagent's sensitivity is 0.22 µg/ml

#### PRO-ZONE EFFECT DETECTED

No pro-zone effect was observed with the high concentration dimer sample up to 90 µg/mL.

#### LINEARITY

The linearity study of the DIMER ONE PLUS product was carried out through the dilution of a concentrated purified D-dimer protein solution. The product's linearity is defined by the highest concentration that the product can measure linearly (through linear regression). The results showed that the product is capable of measuring linearly (through linear regression analysis) up to 70 µg/mL. However, considering that the concentrations of the calibrators, which generate the calibration curve, will be lower than 70 µg/mL, we inform that the product is linear up to the highest point of the calibration curve (highest concentration calibrator). For samples with higher values, it is recommended to dilute the sample with 0.85% Sodium Chloride, repeat the measurement, and multiply the obtained result by the dilution factor. ...

#### CONSUMER INFORMATION:

Biocontinental guarantees the quality of the product, provided that all conditions are respected, and the storage and transport care, specified in the instruction for use, until the expiration date indicated on the packaging.  
If further information or guidance is required, the customer should contact Biocontinental Biological Products.

#### Technical Assistance and request for Instruction of Use in printed format:

Phone: +55 11 40231749 | E-mail: [contact@biocontinental.com](mailto:contact@biocontinental.com)

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Brazilian Industry

**ANVISA 801323900013**

Attention user: check if the version obtained from the Instruction for Use not printed corresponds to the product purchased.

Version: NOV/2025

## UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED